



alō

MEDICAL RECLINER/TRANSPORTER



Shown in treatment height with back style A (610), two tone vanity panel, USB charging port, heat + massage, and fold-away table

Curated for the highest quality of care and aesthetic, the Alō combines form + function to create a stunning selection for any healthcare facility. Completely customizable, the Alō offers premium selections as standard features, allowing maximum flexibility to design the look you need without sacrificing acute functionality. Powered recline, heat + massage, and USB charging port provide modern conveniences and comfort for patients during treatment. Available in two working heights (*treatment* and *recovery*) to accommodate any patient seating need.

STANDARD FEATURES

Manual or powered recline
Swing-away arms
Trendelenburg
Vanity panel selection †
Side table (choice of 2 table styles)
Retractable foot tray
Push handle
Wall-hugger mechanism (infinite recline positions)
Quick-release seat
4" dual thermoplastic casters

† Not available with Tuxedo style



Swing-away arms
open to 180°



Manual release lever or
hand-held power remote



Removable seat for
easy disinfecting

1 ► CHOOSE SEAT HEIGHT



A

STRAIGHT

TREATMENT HEIGHT

RECOVERY HEIGHT

2 ► SELECT BACK STYLE



B

WINGBACK

3 ► VANITY PANEL DESIGNS

Clean and simple
Embossed* (diamond or wave pattern)
Two-tone second material (Grades 1-4 included)


*Available in standard coverings Pearl, Milieu, Dolphin, and Stone.
For more embossed covering selections, contact your local
Champion representative.



Clean and simple



Embossed*



Two-tone

4 ► SIDE TABLE CHOICES

Fold-away or pivot table (available in black only)
--



Fold-away table



Pivot table

STANDARD COVERINGS


 All standard covering selections are Greenguard certified

‡ Iridescent sheen




SAND DOLLAR
C8662

POLYURETHANE




VERDITE
C0592

POLYURETHANE




LUCID
C0460

POLYURETHANE



BALLAD
C0383

POLYURETHANE



DOLPHIN*
MOC249

VINYL



PEARL* ‡
MOC216

VINYL




STONE*
MOC227

VINYL




MILIEU* ‡
MOC210

VINYL




SAGE
SL19

POLY/SILICONE




RAIN
SL30

POLY/SILICONE




GRAPHITE
SL21

POLY/SILICONE



BEACH
SL27

POLY/SILICONE



OCEAN
SL26

POLY/SILICONE

5 ► ADDED CUSTOMIZATION

BODY OPTIONS ◇

Heat + massage	USB charging port (right side only)
Heat only	Articulating tablet arm
Black arm caps	Central-locking casters
Dual-hand control (power models only)	Second side table
Articulating TV/monitor mount	IV pole (2 prong, 4 prong)

◇ Additional charges apply.



USB charging port



Central-locking casters



Black arm caps



Heat + massage



Tablet arm

TUXEDO STYLE ◇

Second material accent on seat back wings, foot rest, and rear cover
Black arm caps
Clean and simple vanity panel
Only available on back style B (wingback)




COVERING UPGRADES ◇

Silicone and Polyurethane
PVC-free material
Graded-in material selections
Customer's own material (COM)

Contact your local Champion rep for more information.



Silica



Patterns




COM

recovery height


SPACE REQUIREMENTS

STANDARD
MODEL


TOP VIEW




TOP VIEW WITH SWING-AWAY ARMS



FRONT VIEW



SIDE VIEW



SPECIFICATIONS

A 610L | 61PL B 630L | 63PL

Powered recline		✓	✓
Manual recline	✓		✓
Weight capacity	300 lbs		300 lbs
Back style	Straight		Wingback
Back height (from floor/from seat)	45"/ 26.75"		48.5"/ 30.5"
Overall width*	32.5"		32.5"
Seat width	23"		23"
Seat height/depth	19.875"/ 20.5"		19.875"/ 20.5"
Arm height (from floor, highest point)**	29"		29"
Overall upright depth (with handle)	38"		38"
Dual thermoplastic casters (diameter)	4"		4"

*Measured with two black fold-away tables.

**Add 1" for model without arm caps

XL MODELS

SPECIFICATIONS

A 620L | 62PL B 640L | 64PL

Powered recline		✓	✓
Manual recline	✓		✓
Weight capacity	500 lbs		500 lbs
Back style	Straight		Wingback
Back height (from floor/from seat)	45.5"/ 27.25"		49"/ 31"
Overall width*	36.5"		36.5"
Seat width	27"		27"
Seat height/depth	19.875"/ 20.5"		19.875"/ 20.5"
Arm height (from floor, highest point)**	29"		29"
Overall upright depth (with handle)	38"		38"
Dual thermoplastic casters (diameter)	4"		4"

*Measured with two black fold-away tables.

**Add 1" for model without arm caps

SPACE REQUIREMENTS

A 620L | 62PL B 640L | 64PL

TOP VIEW

Depth/width	38"/36.5"	38"/36.5"
With swing-away arms	94"	94"

FRONT VIEW

Height	45.5"	49"
Tables down width	36.5"	36.5"
Tables up width	56.5"	56.5"

SIDE VIEW


Upright depth	38"	38"
Trendelenburg depth	75"	78"

treatment height


SPACE REQUIREMENTS

STANDARD
MODEL


TOP VIEW




TOP VIEW WITH SWING-AWAY ARMS



FRONT VIEW



SIDE VIEW



SPECIFICATIONS

A 610 | 61P

B 63P

Powered recline	✓	✓
Manual recline	✓	
Weight capacity	300 lbs	300 lbs
Back style	Straight	Wingback
Back height (from floor/from seat)	47"/ 29"	49"/ 31"
Overall width*	32.5"	32.5"
Seat width	23"	23"
Seat height/depth	24.5"/ 20.5"	24.5"/ 20.5"
Arm height (from floor, highest point)**	29"	29"
Overall upright depth (with handle)	43"	44"
Dual thermoplastic casters (diameter)	4"	4"

*Measured with two black fold-away tables.

**Add 1" for model without arm caps

XL MODELS

SPECIFICATIONS

A 620 | 62P

B 64P

Powered recline	✓	✓
Manual recline	✓	
Weight capacity	500 lbs	500 lbs
Back style	Straight	Wingback
Back height (from floor/from seat)	48"/ 30"	49"/ 31"
Overall width*	36.5"	36.5"
Seat width	27"	27"
Seat height/depth	24.5"/ 20.5"	24.5"/ 20.5"
Arm height (from floor, highest point)**	29"	29"
Overall upright depth (with handle)	43"	44"
Dual thermoplastic casters (diameter)	4"	4"

*Measured with two black fold-away tables.

**Add 1" for model without arm caps

SPACE REQUIREMENTS

A 620 | 62P

B 64P

TOP VIEW

Depth/width	43"/36.5"	44"/36.5"
With swing-away arms	94"	94"

FRONT VIEW

Height	48"	49"
Tables down width	36.5"	36.5"
Tables up width	56.5"	56.5"

SIDE VIEW

Upright depth	43"	44"
Trendelenburg depth	76"	77"

CURATED DESIGN MIX

Need inspiration for further customization of your Alō? Here are some combinations of standard and graded-in coverings, as well as some accessories to provide inspiration for your next design.



Large circles indicate the covering material for the body of the chair, and smaller circles indicate accent materials for two-tone vanity panels, tuxedo styling, or accessories. All of the materials listed, and many other selections, can be found on the Champion website.

championchair.com/fabrics-finishes

champion®

Information contained in this literature is based on the latest product data available at the time of printing. Swatches are enlarged to show detail. Photography, artwork, text and specifications are subject to change. Photography may show optional items. Due to photographic and printing processes, product colors in this brochure are not necessarily accurate. Champion Manufacturing Inc. reserves the right to make changes to colors, materials, prices, standards, options, equipment and any specification at any time without notice or obligation. Contact Champion direct for additional information and the latest specifications.