




Increased padding for added softness and comfort during treatment, subtle curves, and swing-away arms for infection control make the Ascent recliner ideal for both patient and caregiver.

It has all the durable engineering and strength of our classic recliner design with an upscale look.



All steel, powder-coated frame for durability.
Z-spring construction for greater comfort.

STANDARD FEATURES

- Swing-away arms
- Gray fold-away table (one)
- Retractable foot tray
- Push handle
- Wall-hugger mechanism (infinite recline positions)
- Quick-release seat
- 10.2 cm dual thermoplastic casters



Swing-away arms
open to 180°



Quick-release seat




Wall-hugger
recline mechanism

STANDARD VINYL SELECTION


Printed color swatches may not be an exact match to actual fabrics and are enlarged to show detail




NATURAL
T09




AQUAMARINE
T19




BONNIE BLUE
T21




ICED MINT
T26




TEAROSE
T43



COLONIAL BLUE
T45




GRAY
T56




CRANBERRY
T96

SPACE REQUIREMENTS STANDARD MODEL


TOP VIEW




TOP VIEW WITH SWING-AWAY ARMS



FRONT VIEW



SIDE VIEW



BODY OPTIONS

Trendelenburg	USB charging port
Additional gray fold-away table	Heat + massage
Arm caps	Heat only
Central-locking casters	Dual-hand control (power models only)
Cupholder on the fold-away table(s)	

♦ Additional charges apply.



USB charging port




Arm caps



Heat + massage

COVERING UPGRADES

Graded vinyl selections
PVC-free material
Customer's own material (COM)
Silicone and Polyurethane



Graded vinyl selections

RECOMMENDED ACCESSORIES

Articulating tablet arm
Cushion wedge
Head pillow/lumbar support

For a full list of accessories, view our Accessories brochure, the Champion website, or ask your Champion representative for more information.



Articulating tablet arm

SPECIFICATIONS

	STANDARD MODEL NUMBER 65P 650	XL MODEL NUMBER 66P 660
Power recline	✓	✓
Manual recline	✓	✓
Weight capacity	136 kg	227 kg
Back height (from floor/seat)	119.4 cm/73.7 cm	124.5 cm/78.7 cm
Overall width (two gray tables down/up)*	82.6 cm/132.1 cm	92.7 cm/142.2 cm
Maximum width (arms open at widest)	228.6 cm	238.8 cm
Minimum width (arms open at 180°)	82.6 cm	92.7 cm
Seat width	58.4 cm	68.6 cm
Seat height/depth	62.2 cm/52.1 cm	62.2 cm/52.1 cm
Arm height (from floor)	73.7 cm	73.7 cm
Overall upright depth/Trendelenburg depth	109.2 cm/190.5 cm	111.8 cm/195.6 cm
Dual thermoplastic casters (diameter)	10.2 cm	10.2 cm

*Measured with both arms closed.

champion.

Information contained in this literature is based on the latest product data available at the time of printing. Swatches are enlarged to show detail. Photography, artwork, text and specifications are subject to change. Photography may show optional items. Due to photographic and printing processes, product colors in this brochure are not necessarily accurate. Champion Manufacturing Inc. reserves the right to make changes to colors, materials, prices, standards, options, equipment and any specification at any time without notice or obligation. Contact Champion direct for additional information and the latest specifications.