



650 Ford Street
Colorado Springs, CO 80915
Toll: (800) 433-4064 * Ph: (719) 637-0300
Fax: (719) 597-8273 * www.harloff.com

HARLOFF OptimAI™ LINE ALL-ALLUMINUM MEDICATION BOX CART #AL210BOX

Specifications

210 Med Box Capacity (210 left column, 30 boxes x 7 drawers = 210 total)



AL210BOX KEY FEATURES AND BENEFITS:

- Spacious Medication Delivery System with 7 med box drawers, 30 boxes per drawers for 210 med box capacity
- The top of the cart is a smooth and sturdy surface for carrying regulatory documents and pamphlets
- Pull out writing shelf for additional writing space for logging patient's charts
- Easy to push All Aluminum Construction, with Integrated push handles for easy transporting
- Maximum Security, built with BEST® Brand Locks on both cabinet and narcotics drawer
- Ball-bearing, full extension, self-closing drawer slides
- The "slam action" lock on the Narcotics drawer allows the drawer to be closed without using the key. An internal drawer cover keeps the drawer contents secure, even if all other drawers are removed.
- Full wrap around vinyl bumper, protecting the cart
- 5" Sealed ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional, easy to maneuver over all floor types; tile, carpet, linoleum, wood and concrete
- Offered in a "Specialty Package," which includes drawer tray with dividers, for top miscellaneous drawer, waste container and locking sharps container, # **AL210BOXSP**
- Hard-baked powder coat finish

Left Column Drawers, 20.25" W x 17" D:

- 1 ea – 4" Top Multipurpose Storage Drawer with internal divider set
- 6 ea – 3.25" Med Box Drawers
- 1 ea – 9.75" Multipurpose Drawer with 2 adjustable dividers

Twelve available End Panel/Drawer Post color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Dark Blue (-DB), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammertone Dark Blue (-TXDB), Hammertone Gray (-TXGY) and Hammertone Green (-TXGN)

Right Column Drawers, 7.62" W x 17" D:

- 1 ea – 4" Top Multipurpose Storage Drawer
- 1 ea – 9.75" Locking Narcotics Drawer
- 2 ea – 9.75" Multipurpose Storage Drawer

Assembled Cart Dimensions:

- 45.75" H x 36" W x 23.5" D (116.2 cm x 91.4 cm x 59.7 cm)
- Cart Net weight: 157 lbs (71.4 kg)
- Approx. Shipping Weight: 197 lbs (89.6 kg)