

Standard Line Med-Bin Cart with 42 3" Bins Standard Package #SL42BIN3

Specifications



Construction:

- Steel Body, Aluminum drawers

Features:

- Dual Column Medication Delivery Cart with 42 3" wide Medication Bins
- Bins designed for 7-day cassettes
- Dual Locking System complies with regulatory security standards
- BEST® Brand Lock on cabinet; Illinois Lock on narcotics drawer
- The "slam action" lock on the Narcotics drawer allows the drawer to be closed without using the key.
- An internal drawer cover keeps the narcotics drawer contents secure, even if all other drawers are removed.
- Ball-bearing, full extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles
- Full wrap around vinyl bumper, protecting the cart
- 5" Sealed ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional, easy to maneuver over all floor types; tile, carpet, linoleum, wood and concrete
- Hard-baked powder coat finish
- Offered in a **Specialty Package** which includes pullout side shelf, drawer tray with dividers, for top miscellaneous drawer, waste container and locking sharps container - see **SL42BIN3SP**

Left Column Drawers:

- 1 ea – 4" Top Multipurpose Storage Drawer with internal divider set
- 7 ea – 3.25" Med Bin Drawers including 42 3" bins for 7 day cassettes
- 6 ea – 6.5" Multipurpose Storage Drawer

Thirteen available color choices for side panels:

- Hammer Beige (HB), Gloss Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Navy (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)
- Standard Drawer color is Hammer Beige

Right Column Drawers:

- 1 ea – 4" Top Multipurpose Storage Drawer
- 1 ea – 9.75" Locking Narcotics Drawer
- 2 ea – 9.75" Multipurpose Storage Drawer

Assembled Cart Dimensions:

- 45.75" H x 43" W x 23.5" D (116.2 cm x 109.2 cm x 59.7 cm)
- Cart net weight: 210 lbs (kg)
- Approx shipping weight: 275 lbs (125 kg)