

Standard Line/ Model: SL600EPC

Punch Card Medication Cart, key Lock, Standard Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Painted steel cabinet and drawers

Features:

- Dual column medication cart with single width narcotics box with 4 keys
- Basic electronic pushbutton lock, key lock override with 4 keys
- Factory set user and supervisor code
- Lock can be programmed in field
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with intergrated push handles
- Corner bumpers, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish

Left Column Drawer Configuration:

- (1) – 4" multipurpose storage drawer
- (3) - 9.75" punch card drawers

Right Column Drawer Configuration:

- (1) – 4" multipurpose storage drawer
- (1) – 9.75" locking narcotics drawer
- (2) – 9.75" multipurpose storage drawer

Overall Cart Dimensions:

- 45.75" H x 43" W x 23.5" D (116.2 cm H x 109.2 cm W x 59.7 cm D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 195 lbs (88.6 kg)
- Approx. shipping weight: 235 lbs (106.8 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D (144.8 cm H x 119.4 cm W x 68.6 cm D)